



T006

TL Port Configurations

TL's			
TYPE	PORT1	PORT2	
LECTRON TL Corr. D494147 Pulsator 1 x 3 core cable (3mts)	N C	N C	C Pack Corr. D255654MS
MR TL Corr. D494167 Milk Meter 1 x 3 core	N C	N O	C & O Pack Corr. D255655MS
ISOLATOR TL Corr. D494170 ACR, 2pts 2 open cable entries	N C	N C	C Pack Corr. D255654MS
ACR TL Corr. D494177 ACR, 1pt 2 open cable entries	N C	N C	C & O Pack Corr. D255655MS
DVTL Corr. D494175 Arm/Servo 2 open cable entries	N O	N O	O Pack Corr. D255656MS
N C = Normally Closed (Plastic Spindle) N O = Normally Open (Brass Spindle)			